

PHYSIOLASER 500

14 programs

9 Health&Rehabilitation - 5 Beauty&Fitness

Physiolaser 500 accelerates the inflamed tissue recovery, regenerates cells and improves microcirculation.

Physiolaser 500 is especially effective in all the situations where the muscle/articular tissue is inflamed. Inflammation is usually due to a reduced microcirculation. In return the circulation reduction reduces the blood quantity arriving to the cells, causing an ischemic damage.

Reducing the ischemic time, **Physiolaser 500** reduces the pathology duration and the relative pain.

Physiolaser 500 accelerates the healing and the pain elimination both used alone and toghether with other therapeutic modalities. It allows in this way to reduce the use of analgesics and anti-inflammatories.

Physiolaser 500 enhances the natural healing body process without pain.

- 9 Health&Rehab programs
- 5 Beauty&Fitness

Programmable: in multi-spot mode and scansion mode

The main therapeutic effects are:

- Anti-inflammatory effect
- Anti-edemigenous effect
- Biostimulating effect
- Analgesic effect

Physiolaser 500 enhances the natural healing body process.



TECHNICAL FEATURES

Medical Device of II b class

Laser of IV Class

Power (CEI EN 60825-1): 500 mWatt

Beam divergence:13° Laser type: Ga Al As Diode

Beam length: 808

Spot size: 3,5 mm and 11.5mm Operation mode: continuous and

pulsed

Programmable: in multi-spot mode

and scanning mode

In compliance with the directive

93/42/CEE-2007/47/CEE

EQUIPMENT

- 1 Carrying case
- 1 Physiolaser 500 unit
- 1 Handpiece Laser
- 1 pair of glasses for the patient
- 1 pair of glasses for the operator
- 1 Operating manual
- 1 power supply





PHYSIOLASER 500

Program List

Health

Arthrosis

Ankle sprain

Knee pain

Ankle sprain edema

Articular shoulder pain

Epicondilytis

Chronic back pain

Neck pain

Tendinopathies

Beauty

Breast stretch marks

Thighs, hips and gluteus edematous cellulite

Body fat and spread cellulite

Stretch marks

Scares